

CHS Ref. No.: AGS-00101
**METHODS AND DEVICES FOR IDENTIFYING RELATED IONS FROM
CHROMATOGRAPHIC MASS SPECTRAL DATASETS CONTAINING
OVERLAPPING COMPONENTS**

5 **RELATED APPLICATIONS**

This application is a continuation-in-part of Attorney Docket No. 10020515-1
(2003309-0034), U.S. Patent Application No. 10/388,088, filed March 13, 2007, ^{US 20040180446}entitled "
Methods and Devices for Identifying Biopolymers Using Mass Spectroscopy", Dean R.
Thompson and Steven M. Fischer, which is incorporated herein by reference in its entirety.

10

BACKGROUND

Technical Field

This application is related to mass spectral analysis, and more particularly to
15 processing mass spectra generated by mass spectral analysis.

Description of Related Art

Mass spectroscopy is a powerful analytical tool that may be used in identifying
unknown compounds as well as their quantities. Mass spectroscopy may also be useful, for
20 example, in elucidating the structure and chemical properties of molecules, and may be used
in connection with organic as well as inorganic substances. The identification of proteins and
other molecules in a complex mixture derived from biological sources may be performed
using mass spectroscopy. A variety of different techniques have been developed for use with
the identification of molecules, such as proteins.

25

filtered data may be grouped according to the group vector which results from performing the steps of flowchart 700.

5 The example display 1400 of Figure 17 represents the reordered m/z rows such that
m/z rows in the same group are adjacent. After selecting relevant scan(s) for each group, the
corresponding intensities for the selected scans may be obtained from the filtered data set to
produce a resulting spectra. In one embodiment as described herein, the scans may be
selected by finding the scan or time at which each group maximizes the correlation value by
adding the rows of the data matrix for each group and selecting the scan with the maximum
10 intensity value.

The foregoing processing techniques described herein, for example, in connection
with flowchart 400, may not be used in instances where there are two or more molecules that
elute at the same time and also have the same elution profile. In this instance, the foregoing
15 processing steps are not able to identify the different peptides and properly pair parent (U
spectra) with fragments (F spectra), and another processing technique may be used, for
example, as described in Attorney Docket No. 10020515-1 (2003309-0034), AGS-00101
US 20040180446
U.S. Patent Application No. 10/388,088, filed March 13, 2003, entitled "Methods and
Devices for Identifying Biopolymers Using Mass Spectroscopy", hereinafter referred to as
20 "the Thompson and Fischer disclosure". The processing steps of Thompson and Fischer may
be performed on the results produced by processing steps described herein to resolve the
parent-fragment pairings in instances where two or more molecules elute at the same time.
The Thompson and Fischer disclosure describes a method for gathering structural
information for biopolymers in a sample by running the mass spectrometer in the alternating
25 scan mode, as described elsewhere herein, with alternating U and F spectra. Alternating scan